


Multifunctional handset

The handset with built-in microphone offers the patient the possibility to make an alarm call, control lights and sunblinds, television and music.

Remark :

The cable of the handset is not intended for the patient to pull himself up. Avoid therefore connecting the cable above the bed and have the pear hanging. Inform the patient that the push pear is not intended for this purpose.



Functioning of the multi-functional handset

The related button is always indicated between [...] next to the function.

CALL POSSIBILITIES

Making a call [1]

1. Push the red call button for urgent calls. A clear red led informs that a call has been made.
2. The nursing staff can start an intercom connection with the patient via the loudspeaker of the handset. The patient can answer via the microphone on the handset.

Making a service call [7]

1. Push the red call button for service calls. A clear red led informs that a call has been made.
2. The nursing staff can start an intercom connection with the patient via the loudspeaker of the handset. The patient can answer via the microphone on the handset.

LIGHT CONTROL POSSIBILITIES

(Dis)activating the room light [2]

(Dis)activating the reading light [3]

Remark : The weak yellow led near this button helps localising the handset in the dark.

SUNBLINDS POSSIBILITIES (if available)

Opening sunblinds

Push the “open sunblinds” button [5]

Closing sunblinds

Push the “close sunblinds” button [6]

AUDIO and SELECTION POSSIBILITIES

Volume control :

For music, phone and television, volume can be controlled via the button with yellow mark (silence) or red mark (louder).

The audio can be turned off via the “mute” button [16]

Channel control :

Selection via the numeric keypad.

Switch between loudspeaker of the room and on the handset [17]

This function allows to switch to the handset loudspeaker in busy rooms (eg. Two beds rooms, multiple beds, family visits).

TELEVISION POSSIBILITIES

(Dis)activating the television

Push the TV on/off button [11]

Switch between television and teletext

Push the information button [13]

MUSIC DISTRIBUTION

(Dis)activating the music

Push the music on/off button [9]

The channel number is indicated on the display